GERHÓ®

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(GB) HEAT CURING ACRYLIC RESIN WITH A HIGH MOLECULAR WEIGHT FOR COMPLETE, PARTIAL END COMBINED REMOVABLE DENTAL PROSTHESES AND RELINING.

DIRECTIONS:

Heat curing acrylic resin liquid and powder with a high molecular weight for total and/or partial removable protheses, distinguished by long plasticity, and for combined protheses and relinings.

CONTRAINDICATIONS:

- Direct contact of non-polymerised material in the intraoral area.
- In case of proven allergy to components.
- All uses not indicated.

PREPARATION: Mixing ratio: three parts of powder and one part of liquid (2:1 by weight). We recommend working at a temperature of 22 ° / 23 ° C. Pour the liquid into a container with a lid. Slowly add the powder. Mix for at least 30 seconds. Leave the covered mixture to rest for 15 minutes. After this phase, you can remove the mass from the container. Mix the product better using your hands. From here on, the resin will remain plastic for 20 minutes.

FLASK PREPARATION: After eliminating any traces of wax and washing thoroughly with boiling water, let the flask cool to a temperature between 20°C and a maximum of 40°C, and insulate the plaster. Next, pack the prothesis mould with an adequate quantity of mixture. The two parts of the flask are then joined together with one or two pressing operations (at least 2 minutes), using a hydraulic press at a pressure of 120/150 kg/cm².

RECOMMENDED POLYMERIZATION: Place the flask in water at room temperature until completely immersed. Heat the water in about 40/45 minutes at 70 ° C, keep this temperature for 30 minutes, then bring the water to a boil and keep it for 30 minutes, then let it cool slowly in the water for another 30 minutes. Then, removing it from the water, allow the muffle to cool to room temperature.

FINISHING AND POLISHING: Finish as usual, with burs, rubber points, brushes and with pumice. Mirror polish with the help of brushes and polishing products.

N.B.: Observing the above procedures will result in a denture with high stability and excellent tensile strength.

Powder: Store in a cool and dry place

Liquid: Store at a maximum temperature of 30°C

The product is designed for professional use in the dental sector. Follow the instructions for use. The manufacturer shall accept no responsibility for damages arising from different or inadequate use. The user must personally check the suitability of the product for the intended use, particularly if these uses are not reported in the instructions for use. Do not use after expiration. Warning: Any serious accident that may occur with the device must be reported by the user to the manufacturer and the competent authority of the Member State where the user and/or patient are established.

LIQUID: Highly flammable liquid and vapour. Causes skin irritation. May cause an allergic skin reaction. May cause respiratory irritation. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Avoid breathing vapours. Wear protective gloves/protective clothing/eye protection/face protection. IF ON SKIN: Wash with plenty of soap and water. Store in a well-ventilated place. Keep cool. Dispose of contents/container in accordance with local/regional/national/international regulations.

Contains methyl methacrylate, ethylene dimethacrylate.

POWDER: Do not breathe dust. IF ON SKIN: Wash with plenty of soap and water. Contains benzoyl peroxide, methyl methacrylate. May produce an allergic reaction.

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GERHÒ S.p.A. - Via L. Negrelli, 4 - 39100 Bolzano (Bz) - ITALY